




# Biolitec Leonardo<sup>®</sup> Duster Treating Urology Stones

---

LEONARDO<sup>®</sup> Duster  
Super pulsed Thulium fiber laser  
for lithotripsy

**NEW**  
biolitec<sup>®</sup> Thulium  
fiber laser  
for urology



- Future of Lithotripsy
- Superior stone treatment
  - Improved device
  - Ease of use
  - Finer dusting
  - Low retroplulsion

Duration: 30 minutes

Biolitec Leonardo Duster


Participants: HCP's

Presenter: Product & Marketing Team

- I. The Problem
- II. The Treatment Options
- III. Leonardo Duster- Thulium fiber laser from biolitec



Creator: Shidlovski | Credit: Getty Images/iStockphoto



# Climate change and kidney health: Kidney stone formation: implications of global warming

Robert J Fakheri<sup>1</sup>, David S Goldfarb

Affiliations + expand  
PMID: 21451456 DOI: 10.1038/ki.2011.76  
Free article

## Role of Climate Change in Kidney Stone Formation

Ryan Spicard

Affiliation

PMID: 37

Abstract

Kidney

stone

formation


problems

are also

Nephrolithiasis is a common disease across the world that is becoming more prevalent. Although the underlying cause for most stones is not known, a body of literature suggests a role of heat and climate as significant risk factors for lithogenesis. Recently, estimates from computer models predicted up to a 10% increase in the prevalence rate in the next half century secondary to the effects of global warming, with a coinciding 25% increase in health-care expenditures. Our aim here is to critically review the medical literature relating stones to ambient temperature. We have categorized the body of evidence by methodology, consisting of comparisons between geographic regions, comparisons over time, and other factors like sunlight exposure and regional variation in diet that share some contribution, it appears that heat does play a role in pathogenesis in certain populations. Notably, the role of heat is much greater in men than in women. We also hypothesize that the role of a significant human migration (from rural areas to warmer, urban locales beginning in the last century and projected to continue) may have a greater impact than global warming on the observed worldwide increasing prevalence of nephrolithiasis. At this time the limited data available cannot substantiate the hypothesis. Further studies to investigate this effect are warranted.

temperature  
at increase the  
view reports current  
es are more  
mid. Kidney stone  
greater impact on men  
the occurrence of kidney  
increase as well.

## The Stone Belt- Urolithiasis & Nephrolithiasis




**Figure 2:** The so-called stone belt (red) extends all the way around the world and is characterized by urinary stone prevalence of 10 to 15%. In this zone the climatic and social conditions are conducive to stone formation. Some stones are associated with poverty, others with affluence. In Europe and the USA, there has been a sharp, almost exclusively affluence-related rise in the occurrence of calcium oxalate and uric acid stones. Climate simulations for the USA indicate that the stone belt will move northwards in the coming two decades (12).

# The Problem

## Nephrolithiasis or Urolithiasis

- Hard mineral and salt deposits
- Either at Kidney, bladder or ureter
- Calcium stone are the most common one. Uric acid, struvite, and cystine are other types



# The Problem

## Prevalence

- 12% of the world population
- Most common condition affecting the urinary system



# Treatment options

\_\_\_ Medication and drink fluids (<0.5cm)



Intervention:

\_\_\_ Shockwave lithotripsy (ESWL)

\_\_\_ Ureteroscopy/RIRS


\_\_\_ Percutaneous nephrolithotomy (PCNL)

\_\_\_ Open surgery



Reference:  
ESWL Picture: <https://zurology.com/extracorporeal-shock-wave-lithotripsy/>  
PCNL Picture: <https://www.urologyhealth.org/urology-a-z/k/kidney-stones>  
RIRS Picture: <https://www.healthline.com/health/ureteroscopy-with-laser-lithotripsy>

# The Solutions



## Ureteroscopy with laser lithotripsy advantages

- \_\_\_ Effective with high success rate
- \_\_\_ Minimal invasive
- \_\_\_ Can work on multiple stones
- \_\_\_ Same laser for soft tissue application

## The Current standard:

### Holmium YAG laser- Current gold standard

- \_\_\_ High power (140W)
- \_\_\_ Wavelength 2100 nm
- \_\_\_ Safety profile
- \_\_\_ All stone types can be disintegrated



# The New Solutions- Leonardo Duster

## Thulium fiber laser- Revolution in Lithotripsy


- \_\_\_ Low energy per pulse (as low as 0.025J)
- \_\_\_ High frequency range (as high as 2000 Hz)
- \_\_\_ Small core fibers (as low as 150 micron)
- \_\_\_ Uniform beam profile
- \_\_\_ Symmetrical and constant energy pulse
- \_\_\_ Wavelength 1940 nm



# How does it work?

## Thulium fiber laser- Revolution in Lithotripsy

Working principle




# The Solution

## Side by side comparison -Thulium fiber laser and Holmium YAG laser


### Wavelength-

More absorption coefficient in water for Thulium fiber laser compared to Holmium YAG laser.



Source. <https://tau.amegroups.org/article/view/28326/html>

Penetration depth four times lower for Thulium Fiber laser when compared with Holmium YAG laser.




Source. <https://tau.amegroups.org/article/view/28326/html>

# The Solution

## Pulse energy


## Frequency

## Pulse duration and pulse profile



Source: <https://tau.amegroups.org/article/view/28326/html>

# The Solutions



Source. <https://tau.amegroups.org/article/view/28326/html>

# The Solutions



## Advantage overview:

- \_\_\_ Low pulse energy and high frequency
- \_\_\_ Low retro-pulsion
- \_\_\_ Small core fibers
- \_\_\_ Compact device with low noise
- \_\_\_ Air cooling
- \_\_\_ No need of additional basket removal

Parameter	Holmium:YAG laser (Lumenis Pulse 120H)	Thulium fiber laser (IPG Medical, Superpulse)
Wavelength	2120 nm	1940 nm
Pulse energy range	0.2–6.0 J	0.025–6.0 J
Pulse duration range	0.05–1 ms	0.05–12 ms
Pulse shape	Dictated by the pumping pulse	Electronically modulated
Maximum pulse frequency	120 Hz	2000 Hz
Maximum average power	120 W	60 W
Lowest proximal laser fiber core diameter	≥ 200 μm	≥ 150 μm
Cooling system	Low-power generators: self-contained water-cooling system with fan High-power generators: vapor-compression refrigeration system	Fan
Resistance to external shocks	Low	High

Source: Thulium fiber laser: the new player for kidney stone treatment? A comparison with Holmium:YAG laser Olivier Traxer and Etienne Xavier Keller

# Leonardo Duster- Thulium fiber laser from biolitec



- Similar features like other Thulium fiber lasers
- Bare fibers from 150 micron to 1000 microns

## Specifications

Model	LEONARDO® Duster
REF	TFL1940nm60W
Wavelength	1940 nm ± 3 nm
Max average power	60 Watt
Max peak power	600 Watt
Pulse energy	0.02– 6 J
Aiming beam	Laser diode, green 525 nm (adjustable) < 1mW, class 3R (brightness adjustable)
Treatment mode	CW / Pulse
Pulse duration	40 µs to 50 ms
Repetition rate	1 to 2500 Hz
Power supply	110 – 240 VAC, 50/60 Hz
Dimensions (H × W × D)	37.8 cm × 57.7 cm × 36.8 cm
Weight	41 kg

All laser sets incl. 3 safety goggles, foot switch, power cord and manual.  
Trolley or carrying case can be supplied upon request.  
Fiber sizes ranging from 150 to 1000 microns are available.

## Treatment parameters

User defined parameters: Pulse energy, Frequency and Pulse duration (Short, Medium and Long)

Stone treatment is mostly in Short pulse.

\_\_\_ Dusting parameters:

Low pulse energy and High frequency (0.05J-0.2J \* 100-400 Hz)<sup>1</sup>

\_\_\_ Fragmentation parameters:

High pulse energy and Low frequency (0.4J-4J \* 14Hz-40 Hz)<sup>1</sup>

1. Glossary of pre-settings given by laser companies: no consensus! Alba Sierra · Mariela Corrales · Adrià Piñero · Merkourios Kolvatzis · Bhaskar Somani · Olivier Traxer